

# Invia<sup>®</sup> White Foam NPWT

## POLYVINYL ALCOHOL (PVA) DRESSING



### **Effective**

- Open-cell PVA foam ensures reliable removal of wound exudate from the wound.

### **Versatile**

- Can be used for multiple wound types and shapes.

### **Reliable**

- High tensile strength for safe removal without disruption and fragmentation especially in tunneling and undermined wounds.

## Providing choices for Negative Pressure Wound Therapy

The Invia White Foam is an open-cell polyvinyl alcohol (PVA) foam with interlinked pores, optimized for the treatment of wounds by means of negative pressure wound therapy (NPWT). The Invia White Foam is a hydrophilic, effectively absorbent foam. Combined with the Invia Negative Pressure Wound Therapy (NPWT) Systems with negative pressure levels between 80mmHg and 150mmHg, Invia White Foam ensures reliable removal of wound exudate\*. The applied negative pressure is distributed evenly on the base of the wound.<sup>1</sup>

### Key features at a glance

- Packaged moistened with sterile water for the ease of use
- Super-absorbent hydrophilic foam properties makes it less adherent to the wound bed
- Moistening of foam facilitates ease of removal without disruption of granulation tissue formation
- Minimal change of in-growth of granulation tissue facilitates an easy and safe removal



\* Refer to the Invia White Foam NPWT Instructions for Use for details regarding indications and contra-indications.

### Technical data

<b>PVA foam</b> <b>Hydrophilic polyvinyl alcohol foam</b>	Pores per cm: 10–14
	Tensile strength: 194–450 kPa
	Sizes: small (10 x 7.5 x 0.9 cm), large (15 x 10 x 0.9 cm)



The Invia White Foam is intended to be used in conjunction with Invia Foam Dressing Kits with FitPad (available as small, medium, large and X-large).

Find out more about Medela NPWT Portfolio at [www.medela.com](http://www.medela.com)

<sup>1</sup> Data on file at Medela AG.



Medical Vacuum Technology  
for Healthcare Professionals

**Please contact us or your local Medela representative for details.**

Medela and Invia are trademarks of Medela Holding AG, Switzerland.

#### Local contact

**Manufactured by:**  
Mondomed® N.V.  
Middenweg 20  
3930 Hamont-Achel  
Belgium  
TEL. +32 11 44 09 90  
[www.mondomed.be](http://www.mondomed.be)

**Distributed by:**  
Medela AG  
Lättichstrasse 4b  
6341 Baar, Switzerland

**United Kingdom**  
Medela UK Ltd.  
Huntsman Drive, Northbank Industrial Park  
Irlam, Manchester M44 5EG  
United Kingdom  
Phone +44 161 776 0400  
Fax +44 161 776 0444  
[info@medela.co.uk](mailto:info@medela.co.uk)  
[www.medela.co.uk](http://www.medela.co.uk)

**International Sales**  
Medela AG  
Lättichstrasse 4b  
6341 Baar  
Switzerland  
Phone +41 41 562 51 51  
Fax +41 41 562 51 00  
[ism@medela.ch](mailto:ism@medela.ch)  
[www.medela.com](http://www.medela.com)